


## Instruction Sheet BiiB Superstructures

Please follow these instructions in order to minimise errors, ensure prompt turnaround of work and to guarantee that optimum results for your casework are achieved.

### Step Instructions

1. Ensure you have a tissue model prepared so that any implant replicas or abutment replicas are exposed for scanning



Tissue Model



Exposed Abutment Replica

2. All BiiB superstructure designs are to be made from pattern resin, acrylic or a light curable product



Pattern Resin

3. Casework must be packed quite well in order to reduce any damage that may occur during transport
4. Job sheet must be fully completed and signed.

**IF YOU HAVE ANY QUERIES, PLEASE DO NOT HESITATE TO CONTACT STONEGLASS IN THE FIRST INSTANCE ON (+61 2) 9764 1036**